

This Quick Guide is a checklist of all of the steps for performing one coupling reaction using the xMAP Antibody Coupling kit (40-50016) and is not intended to replace the kit user manual (89-00002-00-319). Please refer to the user manual for complete instructions, warnings, and precautions.

Reagent Prep:

- Allow kit to equilibrate to room temperature for 20 to 30 minutes

Reagent Calculations:

- Calculate the volumes required for each reagent:

<p>a. Beads for coupling reaction: Desired number of Beads to couple: _____x10⁶ Concentration of stock vial: _____x10⁶ Beads/mL</p> $\text{Volume of stock needed} = \frac{(\# \text{ of beads to couple})}{(\text{conc of stock vial})}$ $\text{Volume of stock needed} = \frac{(\text{_____ beads})}{(\text{_____ beads/mL})}$	<p>Volume of Bead stock needed for the coupling reaction (Step 4):</p> <p>_____ mL</p>
<p>b. Activation Buffer for Activation step: i. 400 µL for reactions of more than 5x10⁶ Beads ii. 480 µL for reactions of 5x10⁶ Beads or less</p>	<p>Volume of Activation Buffer needed for activation (Step 8):</p> <p>_____ µL</p>
<p>c. Sulfo-NHS for Activation step: i. 50 µL for reactions of more than 5x10⁶ Beads ii. 10 µL for reactions of 5x10⁶ Beads or less</p>	<p>Volume of Sulfo-NHS needed for activation (Step 11):</p> <p>_____ µL</p>
<p>d. EDC for Activation step: i. 50 µL for reactions of more than 5x10⁶ Beads ii. 10 µL for reactions of 5x10⁶ Beads or less</p>	<p>Volume of EDC needed for activation (Step 14):</p> <p>_____ µL</p>
<p>e. Antibody for coupling step: Desired number of Beads to couple: _____x10⁶ Desired Antibody concentration: _____µg/1x10⁶ Beads Stock Antibody concentration: _____µg/mL</p> $\text{Volume of Ab needed} = \frac{(\# \text{ of beads to couple})(\text{Desired Ab conc})}{\text{Stock Ab conc}}$ $\text{Volume of Ab needed} = \frac{(\text{_____} \times 10^6 \text{ beads}) \left(\frac{\text{_____} \mu\text{g}}{1 \times 10^6 \text{ beads}} \right)}{\left(\frac{\text{_____} \mu\text{g}}{\text{mL}} \right) \left(\frac{1 \text{ mL}}{1000 \mu\text{L}} \right)}$	<p>Volume of Antibody needed for coupling (Step 21):</p> <p>_____ µL</p>
<p>f. Activation Buffer for coupling step: Total reaction volume: _____µL i. 1000 µL for reactions of more than 5x10⁶ Beads ii. 500 µL for reactions of 5x10⁶ Beads or less</p> <p>Volume of Activation Buffer needed = total reaction volume – volume of Antibody</p> <p>Volume of Activation Buffer needed = _____µL - _____µL</p>	<p>Volume of Activation Buffer needed for coupling (Step 20):</p> <p>_____ µL</p>

Microsphere Wash #1:

3. Resuspend the stock **Beads**
 - a. 1 mL – vortex and sonicate for 10 seconds
 - b. 4 mL – rotate for 15 minutes at 15-30 rpm
4. Dispense the desired volume of **Beads** (from **Step 2a**) to the reaction tube
5. WASH STEP – use one disposable pipette per reaction tube
 - a. Place reaction tube into magnetic separator (1-2 min)
 - b. Remove supernatant with transfer pipette
 - c. Add 500 µL of **Activation Buffer** into reaction tube
 - d. Vortex and sonicate reaction tube for 10 seconds
6. Repeat WASH STEP, for a total of two washes
7. Remove liquid from **Beads** w/ magnetic separator (1-2 min) and disposable pipette

Activation:

8. Add **Activation Buffer** to reaction tube (from **Step 2b**)
9. Vortex and sonicate reaction tube for 10 seconds
10. Vortex **Sulfo-NHS** tube for 10 seconds
11. Add **Sulfo-NHS** to reaction tube (from **Step 2c**)
12. Add 250 µL of **Activation Buffer** to 10 mg vial of **EDC**.
13. Invert **EDC** vial and vortex for 10-12 seconds
14. Add **EDC** solution to reaction tube (from **Step 2d**)
15. Vortex and sonicate reaction tube for 10 seconds

Incubation #1:

16. Shield from light and rotate the reaction tube for 20 ± 2 minutes @ 15-30 rpm
 - a. Incubation start time: _____
 - b. Incubation end time: _____

Microsphere Wash #2:

17. WASH STEP – use one disposable pipette per reaction tube
 - a. Place reaction tube into magnetic separator (1-2 min)
 - b. Remove supernatant with transfer pipette
 - c. Add 500 µL of **Activation Buffer** into reaction tube
 - d. Vortex and sonicate reaction tube for 10 seconds
18. Repeat WASH STEP, for a total of three washes
19. Remove liquid from **Beads** w/ magnetic separator (1-2 min) and disposable pipette

Coupling:

20. Add **Activation Buffer** to the reaction tube (from **Step 2f**)
21. Add **Antibody** to the reaction tube (from **Step 2e**)
22. Vortex the reaction tube for 10 seconds

Incubation #2:

23. Shield from light and rotate the reaction tube for 2 hours ± 5 minutes @ 15-30 rpm
 - a. Rotation start time: _____
 - b. Rotation end time: _____

Microsphere Wash #3:

24. WASH STEP – use one disposable pipette per reaction tube
 - a. Place reaction tube into magnetic separator (1-2 min)
 - b. Remove supernatant with transfer pipette
 - c. Add 500 µL of **Wash Buffer** into reaction tube
 - d. Vortex and sonicate reaction tube for 10 seconds
25. Repeat WASH STEP, for a total of three washes
26. Remove liquid from **Beads** w/ magnetic separator (1-2 min) and disposable pipette
27. Add 1mL of **Wash Buffer** to the reaction tube
28. Vortex and sonicate reaction tube for 10 seconds
29. Protect from light and store at 2-8 °C until needed