

xTAG® CYP2C19 Kit v3 (I046B0427)

Refer to the product package insert MLD-046-KPI-001 Rev A for complete instructions, warnings and precautions.

A. PCR

1. Thaw reagents at appropriate temperatures.
2. Label appropriate number of 0.2 mL thin-wall PCR tubes.
3. Label one 1.5 mL centrifuge tube as "MM" for Master Mix.
4. Prepare PCR Master Mix:

| 2C19v3 Reagent | x 1 Reaction (μL) |
|--|-------------------|
| DNase, RNase-Free Distilled Water | 5.75 |
| xTAG® 10x HS Taq Polymerase Buffer | 3.00 |
| xTAG® CYP2C19 Kit v3 PCR Primer Mix | 3.00 |
| xTAG® Hot Start Taq, keep at -20°C until ready to use (flick/spin) | 0.25 |

5. Mix/Spin PCR Master Mix.
6. Aliquot 12μL of Master Mix to each PCR tube.
7. Add 3μL of sample/control to each tube. Use dH₂O as Neg (Spin/Mix/Spin). [50ng/μL Total Conc. Optimized]
8. Run cycler (store 2-8°C for up to 24 hrs max).

B. Amplicon Treatment

1. Mix/Spin PCR tubes.
2. Prepare Amplicon Treatment Master Mix:

| 2C19v3 Reagent | x 1 Reaction (μL) |
|---|-------------------|
| xTAG® Shrimp Alkaline Phosphatase, (flick/spin) | 1.20 |
| xTAG® Exonuclease I, (flick/spin) | 0.30 |

3. Mix/Spin Master Mix.
4. Add 1.5μL of the Amplicon Treatment Master Mix to each tube containing PCR product (Spin/Mix/Spin).
5. Run cycler (store at 2-8°C for up to 4 hrs).

C. ASPE

1. Thaw reagents at appropriate temperatures.
2. Label appropriate number of 0.2 mL thin-wall PCR tubes.

3. Label one 1.5 mL centrifuge tube "ASPE MM".
4. Prepare ASPE Master Mix:

| 2C19v3 Reagent | x 1 Reaction (μL) |
|--------------------------------------|-------------------|
| DNase, RNase-Free Distilled Water | 7.40 |
| xTAG® 10x HS Taq Polymerase Buffer | 5.00 |
| xTAG® CYP2C19 Kit v3 ASPE Primer Mix | 4.25 |
| xTAG® Hot Start Taq | 0.35 |

5. Mix/Spin ASPE Master Mix.
6. Aliquot 17μL of Master Mix to each 0.2 mL thin-wall PCR tube.
7. Add 3μL treated PCR product to each tube. Mix/Spin.
8. Run cycler (store at 2-8°C for up to 24 hrs max).

D. Bead Hybridization

1. Thaw reagents below at room temperature.
2. Cut/Label Costar plate.
3. Vortex/Sonicate Bead Mix 10 sec. Repeat.
4. Aliquot 20.0 μL of xTAG CYP2C19 Kit v3 Bead Mix to each well.
5. Transfer 1 μL of ASPE product to corresponding wells. Pipet up/down to mix.
6. Cover with Microseal.
7. Run cycler.

E. Reporter Addition

1. Prepare reporter solution (1:75 dilution):

| 2C19v3 Reagent | x 1 Reaction (μL) |
|---|-------------------|
| xTAG® Reporter Buffer | 74.0 |
| xTAG® Streptavidin, R-Phycoerythrin Conjugate G75 | 1.00 |

2. Pipette 75μL of diluted reporter solution into hybridized samples.
3. Pipette up/down to mix.
4. Incubate at RT for 15 minutes in the dark.
5. Run on Luminex® 100/200™ system.

For accurate results, read the package insert and follow the instructions carefully.

For in vitro diagnostic use only.