

ARIES[®]

Flu A/B & RSV Assay



Welcome to the New Way to Work

The ARIES[®] Flu A/B & RSV Assay is a rapid, accurate method for the detection and differentiation of influenza A virus, influenza B virus, and respiratory syncytial virus (RSV) from nasopharyngeal swab (NPS) specimens using ARIES[®] Systems. Key features include:

- **Comprehensive and Flexible:** The inclusion of multiple pathogens with broad strain coverage provides a single test for pediatric and adult populations. The ability to mask particular pathogens also provides flexibility to meet each physician's order.
- **Excellent Clinical Performance:** High clinical performance reduces the need for confirmatory testing that may be associated with rapid antigen tests, thereby providing results you can trust. An integrated sample processing control ensures the assay run is successful from extraction through amplification.
- **Fast Time to Results:** Answers in less than 2 hours with minimal hands-on time allows for rapid results and may lead to more positive patient outcomes.
- **Reduce User Error:** Internal barcode scanning matches samples to cassettes, enabling Position Independent Results no matter where each cassette is placed. Data input errors are also reduced with electronic ordering through bidirectional laboratory information system (LIS) connectivity.

Performance

All 2,479 eligible prospective clinical specimens were tested by an FDA-cleared molecular comparator and the ARIES® Flu A/B & RSV Assay, as shown below:

Target	Positive Percent Agreement (PPA)		95% Confidence Interval (CI)	Negative Percent Agreement (NPA)		95% Confidence Interval (CI)	"No Call" by Reference
Influenza A	299/312	95.8%	93.0% - 97.8%	2131/2165	98.4%	97.8% - 98.9%	2
Influenza B	45/48	93.8%	82.8% - 98.7%	2417/2431	99.4%	99.0% - 99.7%	0
RSV	270/278	97.1%	94.4% - 98.7%	2165/2201	98.4%	97.7% - 98.9%	0

Refer to Package Insert for additional details: Luminex Corporation | ARIES® Flu A/B & RSV Assay (IVD) Kit Package Insert.

Intended Use

The ARIES® Flu A/B & RSV Assay is a polymerase chain reaction (PCR) based qualitative *in vitro* diagnostic test for the direct detection and differentiation of influenza A virus, influenza B virus, and respiratory syncytial virus (RSV) nucleic acid in nasopharyngeal swab (NPS) specimens from patients with signs and symptoms of respiratory tract infection in conjunction with clinical and laboratory findings. The test is intended for use as an aid in the differential diagnosis of influenza A, influenza B, and RSV in humans and is not intended to detect influenza C.


Negative results do not preclude influenza virus or RSV infection and should not be used as the sole basis for diagnosis, treatment or other management decisions. Conversely, positive results do not rule-out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. The use of additional laboratory testing (e.g., bacterial culture, immunofluorescence, x-ray findings) and clinical presentation must be taken into consideration in order to obtain the final diagnosis of respiratory viral infection.

Performance characteristics for influenza A were established during the 2014-2015 and the 2015-2016 influenza seasons when influenza A/H3N2 and A/H1N1 pandemic were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.

If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to state or local health departments for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

The ARIES® Flu A/B & RSV Assay is indicated for use with ARIES® Systems.

Workflow



The operator simply adds specimen to the sample chamber, puts the cassette in the magazine, loads the magazine into the ARIES® System, and the run will start automatically.

Ordering Information

Product Name	Part Number
ARIES® Flu A/B & RSV Assay Kit	50-10020 (24 tests)
ARIES® Flu A/B & RSV Assay Protocol Kit	CN-0335-01 (one time order only)
ARIES® Two Module System IVD Includes: Instrument System Operation Manual Two Magazines Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand	ARIES-M12V1-IVD
ARIES® M1 System IVD Includes: Instrument System Operation Manual One Magazine Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand	ARIES-M6V1-IVD



ARIES® Systems are crafted to increase laboratory efficiency, ensure result accuracy, and fit seamlessly into today's lean laboratory.



support@luminexcorp.com OR **orders@luminexcorp.com**

For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. ARIES' Systems are class 1(I) laser products. Validation of the LIS compatibility must be performed by the end user.

©2016 Luminex Corporation. All rights reserved. Luminex, xTAG, and ARIES are trademarks of Luminex Corporation, registered in the U.S. and other countries.

HEADQUARTERS

UNITED STATES
+1.512.219.8020
info@luminexcorp.com
www.luminexcorp.com

EUROPE
+31.73.800.1900
europa@luminexcorp.com
www.luminexcorp.com

CANADA
+1.416.593.4323
info@luminexcorp.com
www.luminexcorp.com

CHINA
+86.21.802.31150
info@cn.luminexcorp.com
www.luminexcorp.com

JAPAN
+81.3.5545.7440
info@jp.luminexcorp.com
www.luminexcorp.com