

# NxTAG<sup>®</sup> Respiratory Pathogen Panel + SARS-CoV-2 (CE-IVD)

The NxTAG<sup>®</sup> Respiratory Pathogen Panel + SARS-CoV-2 (NxTAG<sup>®</sup> RPP + SARS-CoV-2), developed for use on the NxTAG-Enabled MAGPIX<sup>®</sup> System, is a qualitative test for the detection of nucleic acids from multiple respiratory viruses and bacteria in nasopharyngeal, oropharyngeal, nasal, anterior nasal, and mid-turbinate nasal swabs, nasal aspirates, and nasal wash in UTM<sup>™</sup>, Liquid Amies (ESwab<sup>™</sup>), or equivalent.

## The NxTAG RPP + SARS-CoV-2 Assay offers:

- **Comprehensive Testing:** Detects 23 pathogens, including SARS-CoV-2 (ORF1ab and M gene), in a single tube, enabling the accurate diagnosis and treatment of respiratory illnesses with overlapping symptoms.
- **Scalable Throughput:** Process up to 96 samples in less than 3 hours post-extraction, accommodating variable, day-to-day testing demand.
- **Minimal Hands-On Time:** Pre-plated, lyophilized reagents facilitate a simple workflow with just one pipetting step, ensuring an easy fit in any lab's daily routine.

## Targets

Viral Targets		Bacterial Targets
Adenovirus	Influenza A 2009 H1N1	<i>Chlamydomphila pneumoniae</i>
Coronavirus 229E	Influenza B	<i>Legionella pneumophila</i>
Coronavirus HKU1	Parainfluenza 1	<i>Mycoplasma pneumoniae</i>
Coronavirus NL63	Parainfluenza 2	
Coronavirus OC43	Parainfluenza 3	
Human Bocavirus	Parainfluenza 4	
Human Metapneumovirus	Respiratory Syncytial Virus A	
Influenza A	Respiratory Syncytial Virus B	
Influenza A subtype H1	Rhinovirus/Enterovirus	
Influenza A subtype H3	SARS-CoV-2	

## Performance

The formulation of the NxTAG RPP + SARS-CoV-2 Assay is identical to NxTAG RPP, with the exception of the additional reagents required for the detection of SARS-CoV-2. No changes have been made to the existing NxTAG RPP reagents, reaction conditions, workflow, or software thresholds; therefore, the performance characteristics of NxTAG RPP are still applicable to NxTAG RPP + SARS-CoV-2. The limit of detection (LoD), analytical reactivity, and specificity (including cross-reactivity) of the panel targets were not impacted by the addition of SARS-CoV-2 to the NxTAG RPP assay.

## Limit of Detection (LoD) of SARS-CoV-2 Tested with the NxTAG Respiratory Pathogen Panel + SARS-CoV-2

The LoD for SARS-CoV-2 in the NxTAG RPP + SARS-CoV-2 Assay was assessed by testing a serial dilution of heat-inactivated SARS-CoV-2 culture fluid (ATCC VR-1986HK, heat-inactivated virus) in pooled negative nasopharyngeal specimens (negative clinical matrix). The LoD titer for SARS-CoV-2 was defined as the lowest concentration at which  $\geq 95\%$  ( $\geq 19/20$ ) of the samples tested generated positive calls. The LoD of the SARS-CoV-2 target in the NxTAG RPP + SARS-CoV-2 Assay is 500 copies/mL.

Based on in silico analysis of each oligo sequence to its binding region in each SARS-CoV-2 sequence, it is predicted that the SARS-CoV-2 sequences available from GISAID as of February 11, 2021—including sequences from the United Kingdom (B.1.1.7), South African (B.1.351 or 20H/501Y.V2), Brazilian (P.1 lineage or 20J/501Y.V3), and Californian (one of five reoccurring mutations that constitute the B.1.429 lineage or CAL20C) variants—are 100% detectable by the NxTAG<sup>®</sup> Respiratory Pathogen Panel + SARS-CoV-2 Assay.


## Clinical Performance of NxTAG<sup>®</sup> RPP + SARS-CoV-2 Assay for SARS-CoV-2 Target

Sample	Number of Samples Tested	Positive	Negative	% Agreement with Reference Method	
Positive	74*	73	1†	PPA	100.0%
Negative	360	0	360	NPA	99.7%
Total	434	73	361		

\*This sample set includes 20 contrived positives in addition to various upper respiratory specimen types (nasopharyngeal swab (NPS), oropharyngeal swab (OP), anterior nasal swab, and nasal aspirate specimens).

†This sample was confirmed positive for SARS-CoV-2 via PCR with bidirectional sequencing.

## NxTAG<sup>®</sup> Workflow (Post-Extraction)



## Ordering Information

Product Name	Part Number
NxTAG <sup>®</sup> Respiratory Pathogen Panel + SARS-CoV-2 (CE-IVD)	I056C0471
NxTAG <sup>®</sup> -Enabled MAGPIX <sup>®</sup> System	MAGPIX-XPON4.1-CEIVD
SYNCT <sup>™</sup> Software	CN-SW47

**Luminex**<sup>®</sup>  
complexity simplified.

**orders@luminexcorp.com or support@luminexcorp.com**

For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence.

©2021 Luminex Corporation. All rights reserved. NxTAG and MAGPIX are trademarks of Luminex Corporation, registered in the US and other countries. SYNCT is a trademark of Luminex Corporation. ESwab is a trademark of COPAN S.P.A.

[luminexcorp.com](http://luminexcorp.com)

HEADQUARTERS

UNITED STATES

+1 512 219 8020

info@luminexcorp.com

EUROPE

+31 73 800 1900

europa@luminexcorp.com

CANADA

+1 416 593 4323

info@luminexcorp.com

CHINA

+86 21 8036 9888

info@luminexcorp.com

JAPAN

+81 3 5545 7440

infojp@luminexcorp.com

SS292674