


## NxTAG® Respiratory Pathogen Panel (RPP) RUO

### Syndromic Testing Made Easy

The NxTAG® Respiratory Pathogen Panel detects 22 respiratory pathogens in an easy to use, single well, closed tube system.

- **Minimal Hands-on Time**  
Simplify your workflow with pre-plated reagents
- **Comprehensive Panel**  
Detect a broad range of targets
- **Scalable Throughput**  
Run 1-96 samples at a time, based on your needs
- **Rapid Turnaround Time**  
Decrease your time to results with a quick protocol
- **Customizable**  
Select any number of targets to customize your panel



#### Viral Targets

Influenza A
Influenza A H1
Influenza A H3
Influenza A 2009 H1N1
Influenza B
Respiratory Syncytial Virus A (RSV A)
Respiratory Syncytial Virus B (RSV B)
Parainfluenza 1
Parainfluenza 2
Parainfluenza 3
Parainfluenza 4
Human Bocavirus
Human Metapneumovirus (hMPV)
Human Rhinovirus / Enterovirus
Adenovirus
Coronavirus HKU1
Coronavirus NL63
Coronavirus OC43
Coronavirus 229E


#### Bacterial Targets

<i>Chlamydomphila pneumoniae</i>
<i>Legionella pneumophila</i>
<i>Mycoplasma pneumoniae</i>

# NxTAG<sup>®</sup> Respiratory Pathogen Panel RUO


Save Time and Streamline Your Workflow

Get Reliable Data with Minimal Hands-on Time



Time estimated on 24 samples run simultaneously.

## Simple and Speedy Workflow




### Step 1

Add 1-96 extracted samples to pre-plated test wells.



### Step 2

Integrated multiplex PCR and bead hybridization.



### Step 3

Read on NxTAG<sup>®</sup>-Enabled MAGPIX<sup>®</sup>.

## Ordering Information

Product Name	Part Number
NxTAG <sup>®</sup> Respiratory Pathogen Panel (RUO)	X051C0451 (96 Tests)

**Luminex**<sup>®</sup>  
complexity simplified.

[orders@luminexcorp.com](mailto:orders@luminexcorp.com) or [support@luminexcorp.com](mailto:support@luminexcorp.com)

For Research Use Only. Not for use in diagnostic procedures. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at [support@luminexcorp.com](mailto:support@luminexcorp.com) to obtain the appropriate product information for your country of residence.

©2015-2021 Luminex Corporation. All rights reserved. Luminex, NxTAG, and MAGPIX are trademarks of Luminex Corporation, registered in the U.S. and other countries.

HEADQUARTERS  
UNITED STATES  
+1 512 219 8020  
[info@luminexcorp.com](mailto:info@luminexcorp.com)

EUROPE  
+31 73 800 1900  
[europe@luminexcorp.com](mailto:europe@luminexcorp.com)

CANADA  
+1 416 593 4323  
[info@luminexcorp.com](mailto:info@luminexcorp.com)

CHINA  
+86 21 8036 9888  
[info@luminexcorp.com](mailto:info@luminexcorp.com)

JAPAN  
+81 3 5545 7440  
[info@luminexcorp.com](mailto:info@luminexcorp.com)

5522441